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

VISA Configure Serial Port


Initializes the serial port specified by **VISA resource name** to the specified settings. The VISA class you wire to the **VISA resource name** input determines the polymorphic instance to use. [Example](#)

Use the pull-down menu to select an instance of this VI.

Select an instance

Place on the block diagram Find on the **Functions** palette

VISA Configure Serial Port

TF

Enable Termination Char prepares the serial device to recognize **termination char**. If TRUE (default), the VI_ATTR_ASRL_END_IN attribute is set to recognize the termination character. If FALSE, the VI_ATTR_ASRL_END_IN attribute is set to 0 (None) and the serial device does not recognize the **termination char**.

U8

termination char calls for termination of the read operation. The read operation terminates when the **termination char** is read from the serial device. 0xA is the hex equivalent of a linefeed character (\n). Change the **termination char** to 0xD for message strings that terminate with a carriage return (\r).

U32

timeout sets the timeout value for the write and read operations.

I70

VISA resource name specifies the resource to be opened. This control also specifies the session and class. Refer to [VISA Resource Name Control](#) for more information.

U32

baud rate is the rate of transmission. The default is 9600.

U16


data bits is the number of bits in the incoming data. The value of **data bits** is between five and eight. The default value is 8.

U8

parity specifies the parity used for every frame to be transmitted or received.

0	no parity (default)
1	odd parity

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Untitled 1 Block Diagram *

File Edit View Project Operate Tools Window Help

13pt Application Font

Globals) (All global members)

Tools

PICKING

VISA Configure Serial Port (Instr).vi Front Panel

File Edit View Project Operate Tools Window Help

13pt Application Font

VISA resource name: COM3 (dropdown menu open with LPT1 and Refresh options)

timeout (10sec): 10000

termination char (0xA = '\n' = LF): \A

Enable Termination Char (T): ON

VISA resource name out: % (dropdown menu)

data bits (8): 8

parity (0:none): None

stop bits (10: 1 bit): 1.0

flow control (0:none): None

error in (no error): status [checked] code 40 source

error out: status [checked] code 40 source

New

- Files
- Projects
- Workspaces
- Other Documents

- ATL COM AppWizard
- Cluster Resource Type Wizard
- Custom AppWizard
- Database Project
- DevStudio Add-in Wizard
- Extended Stored Proc Wizard
- ISAPI Extension Wizard
- Makefile
- MFC ActiveX ControlWizard
- MFC AppWizard (dll)
- MFC AppWizard (exe)
- New Database Wizard
- Utility Project
- Win32 Application
- Win32 Console Application
- Win32 Dynamic-Link Library
- Win32 Static Library

Project name:
AnotherSample

Location:
C:\CURS-11\AnotherSample


- Create new workspace
- Add to current workspace
- Dependency of:
[Empty dropdown]

Platforms:
 Win32

OK Cancel

What kind of DLL would you like to create ?

- An empty DLL project.
- A simple DLL project.
- A DLL that exports some symbols.

< Back

Next >

Finish

Cancel

File Edit View Insert Project Build Tools Window Help

(Globals)

(All global members)

DllMain

```
// AnotherSample.cpp : Defines the entry point for the DLL application.
```

```
//
```

```
#include "stdafx.h"
```

```
BOOL APIENTRY DllMain( HANDLE hModule,  
                      DWORD ul_reason_for_call,  
                      LPVOID lpReserved  
                      )
```

```
{  
    return TRUE;  
}
```

```
|
```

ClassView FileView

Build Debug Find in Files 1 Find in Files 2

(Globals)

(All global members)

testFunction

Workspace 'AnotherSample': 1 project(s)

AnotherSample files

- Source Files
 - AnotherSample.cpp
 - StdAfx.cpp
- Header Files
- Resource Files
- ReadMe.txt

```
// AnotherSample.cpp : Defines the entry point for the DLL application.
//
#include "stdafx.h"

BOOL APIENTRY DllMain( HANDLE hModule,
                      DWORD  ul_reason_for_call,
                      LPVOID lpReserved
                      )
{
    return TRUE;
}

__declspec(dllexport)
void testFunction(long a, long b, long* sum, long* prod)
{
    *sum = a + b;
    *prod = a * b;
}
```

ClassView FileView

Sample.vi Block Diagram *

File Edit View Project Operate Tools Window Help

13pt Application Font

Sample.vi Front Panel

File Edit View

String

C:\C...

Stri


C:\C...

a

0

b

0

Tools

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Connectivity

Libraries & Executables

Libraries & Ex...

Call Library F...

Code Interfa...

System Exec.vi

.NET

Input Device ...

ActiveX

Windows Reg...

Functions

Programming

Structures

Numeric

Comparison

File I/O

Synchronization

Measurement I/O

Instrument I/O

Mathematics

Signal Processing

Data Communication

Connectivity

Express

Favorites

Select a VI...

Array

Boolean

Timing

Waveform

Graphics & So...

Cluster & Vari...

String

Dialog & User...

Application C...

Report Gener...

LabVIEW Help


Hide Locate Back Forward Options

Contents Index Search Favorites

- LabVIEW Help
- Finding Example VIs
- Glossary
- LabVIEW 8.0 Features and Changes
- Activating LabVIEW Software
- Using Help
- LabVIEW Documentation Resources
- Quick Tips
- Getting Started
- Fundamentals
- VI and Function Reference
 - Programming VIs and Functions
 - Measurement I/O VIs and Functions
 - Instrument I/O VIs and Functions
 - Mathematics VIs
 - Signal Processing VIs
 - Data Communication VIs and Functions
 - Connectivity VIs and Functions
 - .NET Functions
 - ActiveX Functions
 - Input Device Control VIs
 - Libraries & Executable VIs
 - Call Library Function Node**
 - Code Interface Nodes
 - System Exec
 - Code Interface Nodes
 - Pipes VIs
 - Port I/O VIs
 - Source Control VIs
 - Windows Registry VIs

Call Library Function Node

Calls a DLL or shared library function directly. The Call Library Function Node supports a large number of data types and calling conventions. You can use it to call functions in most standard and custom-made DLLs and shared libraries. If you want to call a DLL that contains [ActiveX](#) objects, use the [Automation Open](#) VI with the [Property Node](#) and the [Invoke Node](#). The Call Library Function Node is expandable and shows data types for the wired inputs and outputs, similar to the [Bundle](#) function. [Details](#) [Example](#)

Place on the block diagram Find on the **Functions** palette

POLY param 1..n are example input parameters of the library function.

POLY return value is an example return value of the library function.

POLY param 1..n output are example output parameters of the library function.

Call Library Function Node Details

The Call Library Function Node consists of pairs of input and output terminals. You can use one or both terminals. If the node does not return **return value**, the top pair of terminals is unused. Each additional pair of terminals corresponds to a parameter in the parameter list of the function you call in top-to-bottom order. You pass a value to the function by wiring to the left terminal of a terminal pair. You read the value of a parameter after the function call by wiring from the right terminal of a terminal pair.

Right-click the node and select **Configure** from the shortcut menu to display the [Call Library Function](#) dialog box, which you can use to specify the library name or path, function name, calling conventions, parameters, and return value for the node. When you click the **OK** button, the node automatically

Tools

(Globals)

Workspace 'AnotherS

Sample.vi Fro


File Edit View

String

C:\C

Stri

C:

- Visible Items
- Help
- Examples
- Description and Tip...
- Set Breakpoint
- Libraries & Executables Palette
- Numeric Palette
- Create
- Replace
- Configure...
- Create .c File...

AnotherSample.dll - 0 error(s), 0 warning(s)

Call Library Function

Library Name or Path

C:\Curs-11\AnotherSample\Debug\AnotherSample.dll

Browse...

Function Name

testFunction

Run in UI Thread

Calling Conventions

C

Parameter

return type

Type

Void

Add a Parameter Before

Add a Parameter After

Delete this Parameter

Function Prototype:

```
void ?testFunction@@YAXJJPAJ0@Z(long arg1, long arg2, long *arg3, long *arg4);
```

OK

Cancel

Help

(Globals)

LL application.

Call Library Function

Library Name or Path: C:\Curs-11\AnotherSample\Debug\AnotherSample.dll Browse...

Function Name: testFunction Run in UI Thread

Calling Conventions: C

Parameter: return type return type

Type: Void arg1 arg2 arg3 arg4

Parameter Before Parameter After Delete this Parameter

Function Prototype:
void ?testFunction@@YAXJJPAJ0@Z(long arg1, long arg2, long *arg3, long *arg4);

OK Cancel Help

a
0
b
0

AnotherSample.dll - 0 error(s), 0 warning(s)

Call Library Function [X]

Library Name or Path:

Function Name:

Calling Conventions:

Parameter:

Type:

Data Type:

Pass:

Function Prototype:

```
void ?testFunction@@YAXJJPAJ0@Z(long arg1, long arg2, long *arg3, long *arg4);
```